

RESEARCH

Open Access

Characterization of black pigment used in 30 BC fresco wall paint using instrumental methods and chemometry

Tania Gatta^{1*}, Luigi Campanella¹, Carlo Coluzza^{2†}, Vania Mambro¹, Paolo Postorino², Mauro Tomassetti¹, Giovanni Visco¹

From CMA4CH 2010: Multivariate Analysis and Chemometry to Cultural Heritage and Environment Taormina, Italy. 26-29 September 2010

Abstract

Background and methods: Several standard powdered black pigments were characterized by means of thermogravimetry TG-DTG and allied techniques. These pigments were used to make standard plaster frescoes at this purpose prepared. The latter ones were subjected to Raman and reflectance analysis. The results obtained, together with TG data, were chemometrically processed and used to identify an analogous standard fresco fabricated by an unknown commercial black pigment, obtaining excellent results.

Results: The same colorimetric and reflectometric techniques, coupled with suitable chemometric techniques, were then successfully used to identify the type of black pigment present in an ancient roman fresco of the Imperial Age (30 B.C.).

Conclusion: TG-DTG resulted useful techniques to authenticate powdered black pigments. Colorimetry and Raman, but also the only colorimetry, were useful to identify an ancient black pigment in situ.

Aim

The aim of the present research was to identify the black pigment used in a roman fresco dated to 30 B.C. using a non destructive method.

Background

The application required comparison with pre-recorded reflectance spectra of standard plaster frescoes ad hoc prepared.

For the preparation of these “standard frescoes” the previous characterization of utilized different standard black pigments was necessary and was performed by means of different instrumental techniques, particularly, thermogravimetry and differential thermoanalysis (TG-DTG).

On the other hand, in order to characterize and identify a standard fresco prepared using an unknown black pigment, thermogravimetry, Raman microspectroscopy and colorimetry were applied, while to identify the in situ Roman fresco, it was possible to apply only colorimetry and reflectance spectroscopy.

Lastly, Multivariate Analysis techniques were used in order to compare all the available instrumental data.

Results and discussion

Thermogravimetry seems to be the most useful technique to characterize the purchased powdered standard black pigments. The TG and DTG curves values obtained are shown in Fig.1, while the TG data and relative activation energy values (E_a) are summarized in Table 1.

More in detail:

(i) Thermogram of Black-Ivory: step at $T < 150$ °C is due to water loss; step at $T \sim 400$ °C is due to the oxidation of the carbonaceous material, deriving from the

* Correspondence: tania.gatta@uniroma1.it

† Contributed equally

¹Chemistry Department, University La Sapienza, P.le A. Moro, 5, 00185, Rome, Italy

Full list of author information is available at the end of the article

Table 1 Thermogravimetric data and activation energy values (E_a) relative to single steps (all values are the mean of three determinations; for the temperature values $RSD\% \leq 0.5\%$; for the mass loss values $RSD\% \leq 0.1\%$; for the E_a values $RSD\% \leq 0.2\%$). (Continued)

Black-Ivory VS	27		250		555		
	45	3.5	400	14.3	560	4.6	74.9
	245		502		690		
			250	107.1	555	211.2	
			500		710		

proteic material combustion; step at $T \approx 640^\circ\text{C}$ is related to the decomposition of calcium carbonate $\text{CaCO}_3 \rightarrow \text{CO}_2\uparrow + \text{CaO}$, included in the hydroxiapatite, that is carbonated hydroxiapatite [1].

(ii) Thermograms of two different Black Carbon samples: very small step at $T < 100^\circ\text{C}$ is due to the water loss; step at $T \approx 600^\circ\text{C}$ is the only important process due to the oxidation of Carbon to CO_2 , but DTG curves of two samples are different, probably due to the different granulometry of carbon particles. The sample from Z supplier seems more homogeneous, in fact the DTG peak of the oxidation process, at around 600°C , is practically symmetric and more regular than the corresponding peak of sample from supplier P. As a conclusion sample Z has probably a granulometry more homogeneous than sample P. Lastly after the step at 900°C no residue is observed in the both cases.

(iii) Thermogram of black vine: step at $T \approx 50^\circ\text{C}$ shows a fair amount loss of water; step at temperature between 250 and 370°C can be attributed to the oxidation of carbonaceous and organic material, not totally burned, probably coming from cellulosic and ligninic material; step at temperature ranging from 600 to 750°C is due to the decomposition of CaCO_3 . Finally, the steps at around 450°C and 490°C are respectively due, to the presence of hydrated silicates included in black vine pigment during the combustion process of vine wood sample, probably contaminated by silica dust and to the carbon oxidation.

(iv) Finally the steps and the behaviours of TG and DTG curves of the unknown considered sample (i.e. the black-ivory-VS) are very similar to those ones of black ivory standard samples. Looking at the Tab.1 we can also observe that activation energy (E_a) data are more useful than temperature and mass loss data to study the authenticity of a powdered black pigment. Really temperature data are strongly influenced by granulometry and by the mass value of the sample. On the contrary E_a results constant when referred to the same process also even if this process occurs in different samples. At this purpose we can observe the case of the oxidation step in two different black-carbon samples, for which E_a is always ≈ 247 kJ/mol, also when the temperatures of


two steps are relatively different. Similar observations can be performed for E_a values of oxidation process of carbon material in black-ivory and in black-vine. Lastly this consideration appears still very true for values of E_a at the decomposition process of small quantities of calcium carbonates included in the two alone pigments (about 210 kJ/mol) in both cases.

Lastly, Multivariate Analysis techniques, i.e. Principal Component Analysis (PCA) and Hierarchic Cluster Analysis (HCA), were used in order to compare all the available thermogravimetric data [2]. The PCA and HCA representation of these TG data after column centering (see data set in "additional file 1, Tab. A"), are shown in the figures 2(a) and 2(b) respectively. It is interesting that the classification as black-ivory, for the unknown black sample before considered only as "probably", was strongly corroborated [3,4].

EDS and FT-IR data in some cases confirmed the thermogravimetric results (see tables 2 and 3); for instance the presence of calcium carbonate in black-ivory and calcium and magnesium carbonate in black-vine; hydroxiapatite in black-ivory and iron oxides and silica traces in black-vine.

Before identifying the black pigments used in the ancient fresco, it was necessary to know if analytical methods and chemometric techniques used in the present research were capable of effectively distinguishing some top-quality purchased black pigments by which the standard frescoes were prepared. To check the validity of the analytical and chemometric procedures used, all the purchased standard pigments and the standard frescoes prepared by these pigments including the black-ivory-VS pigment, supplied as only "probable" and practically still considered by us as an unknown sample, was tested using three analytical techniques: Thermogravimetry, Colorimetry (Table 4) and Raman microspectroscopy (Table 5) [5].

It can be observed as the Raman microspectroscopy technique proves to be a good tool for recording the presence of any black carbon pigment in frescoes. However the wavelength differences of Raman peaks of different black pigments we studied were found to be in practice not always well obvious (Figure 3) [6].



Briefly, colorimetric, Raman and thermogravimetric data, referring respectively to standard frescoes, including that one of black pigment considered unknown, were digitized using modern smoothing techniques (see data set in “additional file 1, Tab. B”) and processed by chemometric software (see PCA and HCA representation in Figs.4(a) and 4(b)). It was thus possible to identify definitively the unknown black pigment, which can be considered certainly as another black-ivory sample (confidence level 80 %).

Lastly to characterize and identify old Roman black pigment (30 B.C.) only colorimetric data (that is the L*, a*, b* CIELab parameters), obtained by measurement in situ using Minolta software was possible to compare with those ones obtained by similar colorimetric measurements on standard frescoes, using chemometric methods (see Table 6).

Table 2 EDS (Energy Dispersion Spectroscopy) data of standard black-ivory and black-vine samples.

Black-Ivory	C ++++	P +
	O ++	Si n
	Ca +	++
Black-Vine	C ++++	Fe ++
	O ++	Ca +
	Si +++	S n
	Mg ++	Mn n

Legend of concentrations:

++++ = very high
 +++ = high
 ++ = medium
 + = low
 n = negligible

CIELab data of standard frescoes and those ones of old roman fresco were elaborated using PCA and HCA analysis (see data set in “additional file 1, Tab.C”). In practice using the chemometric representation of colorimetric data, displayed in Figs. 5(a) and 5(b), we were able to corroborate that the black pigment of the old Roman fresco under test was probably a black vine pigment. However, in this case it was not possible to include in the data matrix also Raman [7] and thermal analysis data of old roman black pigment, so that the identification in this case had a confidence level of not more than 65 %.

On the other hand the Minolta instrument, being portable, could also allow to record in situ the reflectance spectra in the visible region (400-700 nm), therefore it was possible to compare the reflectance spectrum of the black pigment of the old Roman fresco with the spectra of all well characterized standard black pigments in the standard frescoes (see figure 6 in which it is possible to stress as the reflectance spectrum of the black-vine pigment of the standard fresco almost completely is overlapping that one of the black pigment recorded in the Augustus’ small study room). Indeed, the very small differences between the spectra of ancient and modern pigments are justified by the presence of wax (normally used in Roman frescoes to protect from moisture and to increase colour brightness) and by aging processes [8].

Conclusions

In conclusion, thermogravimetry seems to be the most useful technique to characterize and authenticate powdered black pigments, while, by using only the visible

Table 3 IR (Infra-Red) absorption data of standard black-ivory and black-vine samples.

Black-Ivory [cm^{-1}]	Attribution	Black Vine [cm^{-1}]
	M-O (FeO)	420
	M-O (FeO)	451
	M-O (FeO)	462
550	PO_4^{3-}	
582	PO_4^{3-}	
625	CO_3^{2-}	626
724	CO_3^{2-}	725
873	CO_3^{2-}	
1018	CO_3^{2-}	1017
1107	CO_3^{2-}	

reflectance trend and colorimetric parameters, processed by chemometry, we were still able to identify an ancient pigment in situ. Of course in the lucky case in which all the three techniques (colorimetry, Raman and thermogravimetry) can be simultaneously applied, the identification is more simple and sure. The limitation of the present approach is due to the necessity of record standard frescoes spectra as reference.

Experimental

Samples

The fresco we studied was found on the Palatine in Rome after the excavation of "Emperor Octavianus' house", in the room called "Augustus' small study room (Fig.7)". In the present work the research was restricted to black pigment.

It is well known as in Roman frescoes three different types of black pigments were especially used: black-ivory, black-carbon and black-vine [9,10]. At our recognizing purpose we bought all these pigments. To recognize the old Roman black pigment we took colorimetric and reflectance spectrophotometric measures on standard

Table 5 Raman microspectroscopy data of standard frescoes fabricated using standard black pigments and the considered unknown black pigment (i.e. black-ivory-(VS)).

Identification of characteristic Raman peaks		
	[cm^{-1}]	[cm^{-1}]
Black-Ivory	<i>Experimental</i>	<i>Bibliographic reference (1)</i>
	1346 <i>m</i>	1354 <i>s</i>
	1594 <i>s</i>	1609 <i>s</i>
Black-Carbon P	<i>Experimental</i>	<i>Bibliographic reference (1)</i>
	1337 <i>s</i>	1320 <i>s</i>
	1587 <i>s</i>	1595 <i>s</i>
Black-Carbon Z	<i>Experimental</i>	<i>Bibliographic reference (1)</i>
	1338 <i>s</i>	1320 <i>s</i>
	1590 <i>s</i>	1595 <i>s</i>
Black-Vine	<i>Experimental</i>	<i>Bibliographic reference (1)</i>
	-	1350 <i>w</i>
	1580 <i>s</i>	1573 <i>s</i>
Black-Ivory VS	<i>Experimental</i>	<i>Bibliographic reference (1)</i>
	1349 <i>s</i>	-
	1588 <i>s</i>	-

LEGEND: *s* = strong signal, *m* = medium signal, *w* = weak signal

(1) Values reported by University of Florence. Relative acquisitions were obtained with laser line $\lambda=514.5$ nm, lens 100x, time of acquisition 100", resolution 7 cm^{-1} .

We effected our measures with an instrument Lab Ram Infinity by Jobin Yvon Horiba Group, laser He-Ne ($\lambda=632.8$ nm), filter D06, hole 100, grid 1800 strips/mm, resolution 3 cm^{-1} , acquisition 5x60"

frescoes which we prepared in laboratory (Fig.8) using standard pigments, according to the peculiar Roman technique reported in literature [9]. To obtain the final plaster ("tectorio") to be painted we fabricated several layers of different composition (Fig.9). Our samples consisted of three layers: the first one (arriccio) of pozzolana, lake sand, and slaked lime, the second one (intonaco), opportunely spluttered, consisted of lake sand (thinner) and slaked lime, the last one (intonachino) of few millimeters size, of slaked lime and marble's powder; it was polished to improve the reflectivity [11].

Table 4 Colorimetric (CIELab) data of all standard frescoes fabricated using standard black pigments and the considered unknown black pigment (i.e. black-ivory-(VS)).

Identification of CIELab parameters of standard black pigments in a standard frescoes				
Sample	Illuminating	L	a*	b*
Black-Ivory	D 65	29.71±1.34	0.11±0.02	-0.29±0.09
	D 50	27.44±2.57	0.08±0.03	-0.37±0.16
Black-Carbon P	D 65	33.12±1.11	-0.29±0.03	-3.15±0.24
	D 50	32.87±0.54	-0.76±0.03	-3.33±0.21
Black-Carbon Z	D 65	41.28±0.65	-0.44±0.02	-4.15±0.09
	D 50	40.96±0.61	-1.05±0.04	-4.24±0.15
Black-Vine	D 65	29.15±0.61	0.16±0.02	0.38±0.14
	D 50	29.37±0.46	0.21±0.03	0.35±0.12
Black-Ivory VS	D 65	29.41±0.34	0.21±0.04	0.35±0.14
	D 50	27.81±0.46	0.23±0.03	0.44±0.12

Measurements performed by spectrophotometer Minolta CM2600d, standard observatory 2°, MAV mask. Illuminants: D 65 (European standard) and D 50 (American standard).




Figure 3 Raman spectroscopy curves of standard frescoes fabricated using standard black pigments and the considered unknown black pigment sample (i.e. black-ivory-VS).

Standard black pigments were purchased over the counter in Rome and Florence in shops specializing in the sale of restoration materials. The standard black pigments available on the market and used to construct the training set in the chemometric process were: one black-ivory, two black-carbon and one black-vine pigments. To reach our aim we used these samples because they are the only certified samples available on the Italian market.

Lastly, another black pigment, supplied to us (but “only probably”) as black-ivory, was considered as an unknown sample and used for the “Validation Set”, therefore name black-ivory-VS.

Thermogravimetry

We used thermogravimetry (TG, DTG) to characterize standard powdered pigments. The thermogravimetric (TG) curves show the variations of the percent mass of

the sample as a function of increased temperature, therefore the heating generally can determine variations (losses) of the sample mass. The temperatures of these processes are typical for each sample [12]. The DTG curves show the first derivative of the TG curves. These measures were carried out by a thermobalance Mettler TG 10-TA in the following conditions: dynamic air, flow rate 10 ml/min, heating rate 10°C/min from 25 to 900°C, the sensitivity of the balance was 1 µg. It was experimentally determined also the activation energy (E_a) of the main TG steps. E_a was calculated from a single dynamic thermogravimetric measurement using a multiple linear regression (n^{th} order kinetics) method (i.e. the Wyden-Widmann’s method [13]), by means of the following Arrhenius type equation: $d\alpha/dt = k_0 e^{-E_a/RT} (1-\alpha)^n$. In the applied least square method the sum of the squares of the differences between the $d\alpha/dt$ values obtained from the above mentioned equation and values derived from the TG and DTG measurements, respectively taking into account that: $\alpha = \Delta m/\Delta m_{\text{tot}}$ and $d\alpha/dt = (dm/dt)/\Delta m_{\text{tot}}$, attains its minimum values for the selected values of k_0 , n and E_a .

EDS and FT-IR

Using an energy dispersion spectrometer (EDS) as detector we obtained the chemical analysis of standard pigments by scanning electron microscopy, measuring the energy and distribution of the intensity of X photon generated by the electronic beam of the (SEM) electron microscopy, of which the sample under analysis is the target. Atoms having atomic number lower than that one of sodium are not revealed. The IR spectra were obtained by a spectrophotometer FT-IR Perkin-Elmer NIR. At this purpose 2-3 mg of each pigment were dispersed in 100 mg of anhydrous KBr, then obtaining tablets to be analyzed.




Figure 4 (a): PCA representation of the scores of data set values reported in Tab. 6, i.e. CIELab, TG-DTG and Raman data of all the standard black pigment frescoes fabricated in our laboratory and of the considered unknown black pigment sample (i.e. black-ivory-VS). (b): HCA representation of the data set values reported in Tab. 6, i.e. CIELab, TG-DTG and Raman data of all the standard black pigment frescoes fabricated in our laboratory and the considered unknown black pigment sample (i.e. black-ivory-VS).

Table 6 Comparison of colorimetric CIE Lab data of standard fabricated frescoes using standard black pigments and the same data found for the black pigment of the old Roman fresco.


Identification of a black pigment used in the old Roman fresco by means of CIE Lab parameters (values are the mean of three determinations)				
Sample	Illuminating	L	a*	b*
Black-Ivory	D 65	29.71±1.34	0.11±0.02	-0.29±0.09
	D 50	27.44±2.57	0.08±0.03	-0.37±0.16
Black-Carbon P	D 65	33.12±1.11	-0.29±0.03	-3.15±0.24
	D 50	32.87±0.54	-0.76±0.03	-3.33±0.21
Black-Carbon Z	D 65	41.28±0.65	-0.44±0.02	-4.15±0.09
	D 50	40.96±0.61	-1.05±0.04	-4.24±0.15
Black-Vine	D 65	29.15±0.61	0.16±0.02	0.38±0.14
	D 50	29.37±0.46	0.21±0.03	0.35±0.12
Black old Roman fresco	D 65	28.81±0.89	0.42±0.05	2.28±0.17
	D 50	28.54±0.34	0.52±0.04	2.88±0.19


Measurements performed by spectrophotometer Minolta CM2600d, standard observatory 2°, MAV mask. Illuminants: D 65 (European standard) and D 50 (American standard).

Raman microspectroscopy

For Raman microspectroscopy used to characterize our standard frescoes a Lab Ram Infinity (Jobin 1024x256 pixels, Yvon Horiba Group) was used with 25 mW of maximum CCD (Charge Couplet Device) light stimulation, origin of an electric signal. The laser attenuation is given: 1-2% by mirrors, 1-2% by lenses, 50% by beam splitter, 5% by the lens and for a 10^{-6} factor by the notch filter. The laser beam power may be reduced by the regulation of opportune filters which are in order: (starting from the filter with maximum attenuation) D3-D2-D1-D06-D03-without filter (according the signal attenuation

factor due to the selected filter, for example the D03 filter reduces the signal power of a factor 10^3). The high intensity of the laser beam produces a local overheating which can damage the sample modifying its structure. In the case of black pigments (essentially combustive materials) the risk is minimum and it is sufficient to use a filter D03. For the instrument calibration neon was used. For black pigments the only interesting region is the one included in the frequency range of $1200 - 1600 \text{ cm}^{-1}$ (8333.3-6250 nm). We utilized a grid with 1800 strips/mm pointed on the frequency of 1500 cm^{-1} (6666.7 nm), D03 filter, 50x lens, hole 100. The definitive spectrum







was the average of five measures, each one being recorded in 60”.

Visible reflectance spectroscopy and colorimetry

To identify and characterize the black pigments of frescoes we used the colorimetric and reflectometric techniques. At this end colorimetric spectrophotometer Minolta CM2600d allowed us to perform reflectance measurements in the visible region (400-700 nm). It was so possible to determine both the reflectance curve and the values of colorimetric parameters. We measured each commercial pigment on 12 different zones on the standard fresco, obtaining 12 reflectance spectra, of which we calculated the average and the relative standard deviation. It is well known that there are many different ways to define each colour; the more used is the CIE Lab [14], by which each is defined by means of three parameters: L^* , a^* , and b^* . These three parameters are in mathematical connection with Munsell parameter defined as L^* (the luminance), $H^* = \arctan(b^*/a^*)$ (the hue), and $C^* = (a^{*2} + b^{*2})^{1/2}$ (the saturation or chromaticity). In order to define a true colour we used an



“observer of 2°”, that corresponds to cones vision, defined by CIE (Commission Internationale de l’Enclaire) since 1931.

Chemometric methods

For chemometrics the software used was: Lotus spreadsheet 9.8 (IBM/Lotus, Usa); Past 2.14 (free version by Øyvind Hammer, Norway); Datalab 2.99 (light version by H. Lohninger, Austria); Multivariate Analysis (an Excel ad-in by Prof. R. G. Brereton).

Additional material

Additional file 1: Table A: Table_A., *.doc, data set of TG-DTG data after col. cent, “data set of TG-DTG main data of standard powdered black pigments, used to corroborate the identification of black-ivory -VS pigment considered as unknown pigment.” Table B:

Table_B, *.doc, data set of CIELab, TG-DTG and Raman data, "data set of CIELab, TG-DTG and Raman data, used to identify (after column centering) the "black-ivory-VS" considered an unknown sample." Table C: Table_C, *.doc, data set of colorimetric CIELab data, "data set of colorimetric CIELab data used for the identification of the type of black pigment of old Roman fresco."

doi:10.1186/1752-153X-6-S2-S2

Cite this article as: Gatta *et al.*: Characterization of black pigment used in 30 BC fresco wall paint using instrumental methods and chemometry. *Chemistry Central Journal* 2012 **6**(Suppl 2):S2.

Acknowledgements

This work was funded by University of Rome "La Sapienza", Ateneo and University Projects.

This article has been published as part of *Chemistry Central Journal* Volume 6 Supplement 2, 2012: Proceedings of CMA4CH 2010: Application of Multivariate Analysis and Chemometry to Cultural Heritage and Environment. The full contents of the supplement are available online at <http://journal.chemistrycentral.com/supplements/6/S2>.

Author details

¹Chemistry Department, University La Sapienza, P.le A. Moro, 5, 00185, Rome, Italy. ²Physic Department, University La Sapienza, P.le A. Moro, 5, 00185, Rome, Italy.

Authors' contributions

TG collected data and helped draft the manuscript, LC coordinated the research group, CC dealt with colorimetric data, VM performed colorimetric and Raman measurements and fabricated standard frescoes, PP dealt with Raman data, MT performed thermogravimetric measurements and wrote the paper and GV dealt with chemometric processing of data.

Competing interests

The authors declare that they have no competing interests.

Published: 2 May 2012

References

1. Becker A, Epple M, Muller KM, Schmitz I: **A comparative study of clinically well-characterized human atherosclerotic plaques with histological, chemical, and ultrastructural methods.** *J Inorg Biochem* 2004, **98**(12):2032-2038.
2. Wold S: **Chemometrics what do we mean with it and what do we want from it.** *Chemometrics and Intelligent Laboratory Systems* 1995, **30**:109-115.
3. Lavine BK: **Clustering and Classification of Analytical Data.** In *Encyclopedia of Analytical Chemistry*. John Wiley & Sons Ltd, Chichester; Robert A. Meyers 0471-97670-91-20.
4. Gurden SP, Westerhuis JA, Bro R, Smilde AK: **A comparison of multiway regression and scaling methods.** *Chemometrics and Intelligent Laboratory Systems* 2001, **59**:121-136.
5. Vandenaabeele P, Edwards HGM, Moens L: **A Decade of Raman Spectroscopy in Art and Archaeology.** *Chemical Reviews* 2007, **107**(No. 3):675-686.
6. Smith DC, Barbet A: **A Preliminary Raman Microscopic Exploration of Pigments in Wall Paintings in the Roman Tomb Discovered at Kertch, Ukraine, in 1891.** *J. Raman Spectrosc* 1999, **30**:319-324.
7. Villar SEJ, Edwards HGM: **An extensive colour palette in Roman villas in Burgos, Northern Spain: a Raman spectroscopic analysis.** *Anal Bioanal Chem* 2005, **382**:283-289.
8. Ball P: **Colore una biografia.** BUR Saggi 9788817868129 2004.
9. Colombo L: **I colori degli antichi.** Nardini, Roma 9788840440354 1995.
10. Matteini M, Moles A: **La chimica del restauro.** Nardini, Roma 2002.
11. Augusti S: **I colori pompeiani.** De Luca, Roma 1967.
12. Campanella L, Tomassetti M: **Metodi di analisi termica (parte seconda).** *Inquinamento* 1987, **6**:67-69.
13. Wyden H, Widman G: **Kinetic evaluation of TG measurements with multiple regression.** *Angew Chem Thermodyn Thermoanal* 1979, **37**:284-290.
14. Palazzi S: **Colorimetria, La scienza del colore nell'arte e nella tecnica.** Nardini, Firenze 8840440402 2000.

Publish with **ChemistryCentral** and every scientist can read your work free of charge

"Open access provides opportunities to our colleagues in other parts of the globe, by allowing anyone to view the content free of charge."

W. Jeffery Hurst, The Hershey Company.

- available free of charge to the entire scientific community
- peer reviewed and published immediately upon acceptance
- cited in PubMed and archived on PubMed Central
- yours — you keep the copyright

Submit your manuscript here:
<http://www.chemistrycentral.com/manuscript/>


ChemistryCentral